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ADDRESS FOR LETTERS: KIPP DELFT HOLLAND

TELEGRAPHIC ADDRESS: ANKERSMIT DELFT HOLLAND

N.V. INSTRUMENTFABRIEKEN -HANDEL  
**V.H. P. J. KIPP & ZONEN.**  
 DELFT - HOLLAND

---

## Kipp's special lantern for tungsten Arc-lamps.

Philips' glowlampworks at Eindhoven (Holland) are manufacturers of a self-starting tungsten arc-lamp (fig. 1). Two tiny spheres of tungsten, placed closely together, are heated to a high temperature by an arc between the spheres.

A luminous discharge, produced in the bulb when switching on the lamp, ionizes the rare gas with which the bulb is filled and thereby starts the arc.


Fig. 1.

An alternating current of at least 190 V. is necessary to start this lamp, but the tension at the arc is about 30 V. only, the difference of voltage being absorbed by the resistance, which is supplied with the arc-lamp. A special transformer can also be supplied to reduce the loss of energy in the resistance, and to permit the use of the tungsten arc-lamp on A. C. supply of less than 190 V.

The essential quality of the tungsten arc-lamp is its high intrinsic brilliancy, which is about 13000 C. P./sq. inch (21 C. P./mm<sup>2</sup>). This high brilliancy, combined with the absence of need for attendance (necessary when using an ordinary arc-lamp), makes this lamp specially suitable for microphotography, microprojection, dark ground microscopy, galvanometer- and oscillograph-work, etc. To facilitate the use of the tungsten arc-lamp for these purposes, we build a special lantern of tubular form, intended to contain a horizontal type tungsten arc-lamp. The lantern, mounted on a cast-iron foot, can be adjusted in any direction and to various heights. The foot carries a socket for the resistance, which is made in lampshape; 6 ft. of flexible cord are supplied with the lantern, which on account of its small current consumption (1.3 Amp.) can be connected to any electric outlet.

When used with a lens of 30 mm focus, the lantern produces an image of 50 mm in diameter of one of the tungsten spheres at one meter distance. The tungsten-spheres lying, in contradiction of the illustration (fig. 1), both in the optical axis of the lantern, one of them only is effective, the second being shielded from view by the first sphere.

The lamp can be adjusted in the lantern in any direction. When using a lens of 30 mm focus, which is specially recommended for the illumination of mirrors in microscopy, the lamp should be brought very near to the lens.

For galvanometers on the contrary the beam of light ought to be less concentrated for this purpose a lens of 60 mm focus will serve best. The lamp should then be moved to the back of the lantern.

Diaphragms with a circular hole and a slit will be found useful when using this lantern. Both are supplied with the lantern and are to be placed in front of the lens. Two lenses


Fig. 2.

of 60 mm focus are included in the outfit; when using both together the focus is 30 mm.

As stated, a special transformer with resistance and change-over-switch can be supplied with these tungsten arc-lamps.

## Prices.

Special adjustable lantern for tungsten arc-lamps on a cast iron foot provided with resistance, two lenses of 60 mm focus each, two diaphragms, flexible cord and wall-plug, complete with tungsten arc-lamp (for alternating current of 190, 200, 210, 220, 230, 240

250 Volts and higher, the voltage to be stated when ordering) . . . . .	65.—	dutch guilders.
Spare tungsten arc-lamp for alternating current 1,3 Amp. . . . .	7.50	" "
Spare resistance in lamp-form (mention voltage of supply when ordering) . . . . .	3.—	" "

The transformer is built on a wooden base together with 1<sup>o</sup> a socket for the Philips' lampshape-resistance, 2<sup>o</sup> damping resistance for the arc-lamp, 3<sup>o</sup> starting switch, which will always be moved into the running position by a spring, and 4<sup>o</sup> a wall-contact for connecting the lantern. The outfit is provided with 2 Meters flexible cord with plug, and can be connected to outlets of 100 till 140, 190/240 Volts.

Price . . . . . 45.— dutch guilders.

When using this transformer-outfit, the lamp-shape resistance must be screwed into the socket on the wooden base, and a short circuiting plug must be screwed into the socket on the cast-iron plate of the lantern.

---

### OTHER LEAFLETS:

- Galv 24. Moll coil-galvanometers.
- TS-Z 24. Torsion-string galvanometer for direct- and alternating current measurements.  
Zernike coil-galvanometer.
- Acc 24. Accessories and applications for galvanometerwork as: galvanometer-lanterns, scales, shunt-boxes, reflecting prisms, switch-keys, resistances, thermo-elements, alternating current-converter.
- Ther 24. Moll thermopiles for radiation and spectral work.
- Mf 24. Moll registering microphotometers.
- Ex 24. Moll nephelo- and absorptiometer.
- Wolf 24. Tungsten arc-lamp lantern.
- Reg 25. Photographically registering drums.
- WAV 25. Precision alternating current instruments for measuring Volts, Amperes, Watts and Power-factors.
- \*Rapyr 25. Radiation pyrometers.
- \*Bico 25. Gerretsen's Bi-colorimeter for Ph. investigations.
- \*Diaf 25. Noyon's Diapherometer.
- \*Mono 25. Van Cittert double monochromator.
- \*Vibra 25. Moll vibration galvanometer.
- \*Vacuum 25. Moll vacuum thermo-element and vacuum thermo-relais.

\* Will be published within a few months.

ADDRESS FOR LETTERS: KIPP DELFT HOLLAND  
 TELEGRAPHIC ADDRESS: ANKERSMIT DELFT HOLLAND

# N.V. INSTRUMENTFABRIEK EN -HANDEL v.h. P. J. KIPP & ZONEN. DELFT - HOLLAND

## Radiation Pyrometers.

The principle applied in our radiation pyrometers is based on the Stefan Boltzmann's law, being: the total radiation-energy of a black body is proportional to the 4th power of the absolute temperature.


Fig. 1.

So by measuring the total radiation one can find the temperature of the radiating body. This measuring of the radiating energy is accomplished with a Moll-thermopile. The 18 elements of this thermopile consist of thin strips of constantan and manganine of about 3/10 mm wide

and .007 mm thick, soldered together. They are placed behind a circular opening of about 6 mm diameter. Owing to the strips being very thin, their heat capacity is small and the final value of the thermo-electric force is quickly reached.

For scientific work in laboratories etc. we construct a pyrometer type W, as illustrated by fig. 1. This instrument consists of a tube with diaphragms, in one end of which the thermopile is mounted, while the other end is provided with a water-cooled diaphragm, which may be removed if desired.


Fig. 2.

The stand is adjustable and the telescope can rotate horizontally and vertically.

A small window is introduced in the thermopile in order to simplify the sighting operation. To eliminate reflexes the tube is lined with a dull material and diaphragms are provided.

We construct three types of pyrometers for industrial purposes.

**Type A** consists of a tube of fire proof material, built into the furnace by means of a fireclay brick. The pyrometer, which is pointed to the closed end of this tube, is placed on a wall-bracket (fig. 2), fastened to the furnace. A system of lenses concentrates the radiation on to the front of the thermopile.

An accompanying millivoltmeter is provided with a temperature scale reading approx. from 600 to 1450° C. (fig. 3).

For a second range (900° C. and over) an extra diaphragm should be fixed before the lens.

In fig. 2 we have indicated the principal dimensions of the parts which must be fixed to the furnace. In view of the different wall-thickness of furnaces, we supply the pyrometers

with tubes of various lengths. The tubes are made of fire proof material and are able to stand any temperature up to 1500° C.


Fig. 3.

Fig. 4 illustrates the pyrometer **type B**, which differs from the previous model in as much that the thermopile and the system of lenses is mounted in a sheet iron case, fixed on a base with bent border, in which is a hole through which the flange of the fire proof tube can pass. However this instrument can be supplied with one temperature scale only.


Fig. 4.

In the iron case is a small window with a protection wire-gauze; this window may be shut off by a drop-slide, suspended by a lead wire.

This type is intended for use on furnaces subject to overpressure, which causes the flames to escape through the aperture, if the fire proof tube should accidentally break down. Immediately this happens, the lead wire fuses and the slide drops, thereby protecting the pyrometer from damage by overheating.

Both type A and B have lenses which not only concentrate the radiation, but also absorb a part of it and since this absorption varies according to the temperature of the furnace. The application of the Stefan Boltzmann's law does therefore **not** indicate the relation between thermo-force and temperature. For this reason these pyrometers should be specially calibrated in our laboratories. They may also be connected to a recording millivoltmeter, for which we recommend one of the following recorders.

1). A simple clock-work driven recorder of Cambridge Instrument Comp. Ltd. The width of the 48 hour charts is 100 mm and they run with a speed of 12,5 mm (5 inches) per hour. The instrument is mounted in a metal case and is fitted with a temperature scale.

Fig. 5.

2). Fig. 5 illustrates a recording instrument from the same firm, also in metal case with temperature scale. It is also clock-work driven and the 25 hour chart, 100 mm wide, runs with a speed of 12,5 mm per hour. It can be set so that one, two or three independent and differently coloured records can be made simultaneously on one chart. This instrument can also be supplied with two recording mechanisms, in order that 4 to 6 records on one chart of double width may be obtained.

**Type C** (fig. 6) is a portable instrument without lenses, consisting of a metal tube with diaphragms and a thermopile at one end. A millivoltmeter indicates the thermo-force


Fig. 6.

generated. As no lenses are used in this pyrometer, the indication is not affected by the distance to the source of radiation.

We recommend type C for the measuring of temperatures of hardening, annealing and boiler furnaces, also of retorts and such like. It is also suitable for establishing the temperature of molten metals, but in that process the influence of the rather rapidly formed oxyde layer must not be overlooked when calibrating the instrument.

It has proved to be necessary to protect the sensitive part of the thermopiles of types A and B and this is effected by a glass plate. With type C a plate of fluorite is provided, which is a non-selective material to radiation from bodies with a temperature of 800° C. or over.

For some years a great number of our pyrometers have been in continuous use at Philips' Glowlampworks. Dr. Holst, the Manager of the physical laboratories of that firm, after having tested a radiation pyrometer type C, wrote as follows.

„Herewith I give the thermo-force of the radiation pyrometer tested by me at various temperatures, with which I must point out that the diameter of the thermopile area exposed to radiation was 6 mm and that there was a circular diaphragm with an opening of 6 mm at a distance of 22 cm.

$t^{\circ}$	$E \cdot 10^{-6} V.$	$\frac{E \cdot 10^{-6}}{T^4}$
1042	316.1	1.058
1358	756.9	1.068
1485	1017.1	1.065
1558	1202.1	1.068

You see that Stefan Boltzmann's law may be applied, so that calibration in one single point is sufficient. At a distance of about 20 cm one can measure the temperature of a little furnace of 2 cm in diameter with it".

We respectfully request our customers to furnish us with full particulars as to their requirements, when ordering, viz.:

1. Object.
2. What type is required (W, A, B or C).
3. Maximum temperature.

Type A is normally supplied for 600 to 1450° C and type B for 600 to 1100° C, or 900 to 1450° C, or for one specific range as desired, providing there is no greater range than 500° C and the lowest temperature above 600° C. Type C and W can be supplied for any maximum temperature.

4. For type A and B the wall thickness of the furnace.
5. What measuring instrument is required.

We advise the use of the millivoltmeter illustrated in fig. 6 for transportable outfit. The millivoltmeter fig. 3 or one of the recorders should prove more suitable for permanent installations.

The indicators of Cambridge Instrument Comp. Ltd. are enclosed in strong, dust-proof metal cases and are perfectly suitable for continuous workshop use.

For type W the indicator to be used should be selected in accordance with each particular purpose. As this type is generally used to control the temperature of small furnaces in laboratories, this pyrometer has a small opening, so that for low temperatures only a galvanometer will indicate.

**PRICES:**

	Dutch guilders.
Pyrometer type W . . . . .	. . . . .
Pyrometer type A with chamotte tube and fireclay brick . . . . .	. . . . .
Pyrometer type B with base . . . . .	. . . . .
Pyrometer type C . . . . .	. . . . .
Millivoltmeter as illustrated in fig. 3 . . . . .	. . . . .
Millivoltmeter as shown in fig. 6 . . . . .	. . . . .
Recorder as described on page 4, pos. 1) . . . . .	. . . . .
Recorder as described on page 4, pos. 2) fig. 5 . . . . .	. . . . .
For 1 record 167/41184 . . . . .	. . . . .
For 2 records 167/41185 . . . . .	. . . . .
For 3 records 167/41186 . . . . .	. . . . .
For 4 records 167/41195 . . . . .	. . . . .
For 6 records 167/41196 . . . . .	. . . . .
Spare chamotte tube for type A . . . . .	. . . . .
Spare fireclay brick for type A . . . . .	. . . . .

## Thermo-electric measuring instruments.

We manufacture thermo-electric Ampère-, Volt- and Wattmeters as precision instruments, which have the great advantage that they are independent of phase-displacement and frequency. Moreover they can be calibrated with direct current. The wattmeter can specially be made for measuring dielectric losses of cables.

Ask for our leaflet WAV 25.

---

### OTHER LEAFLETS:

- |             |  |
|-------------|--|
| Galv 24.    | Moll coil-galvanometers.   |
| TS-Z 24.    | Torsion-string galvanometer for direct- and alternating current measurements, Zernike coil-galvanometer.   |
| Acc 24.     | Accessories and applications for galvanometerwork as: galvanometer-lanterns, scales, shunt-boxes, reflecting prisms, switch-keys, resistances, thermo-elements, alternating current-converter. |
| Ther 24.    | Moll thermopiles for radiation and spectral work.  |
| Mf 24.      | Moll registering microphotometers.   |
| Ex 24.      | Moll nephelo- and absorptiometer.  |
| Wolf 24.    | Tungsten arc-lamp lantern.   |
| Reg 25.     | Photographically registering drums.  |
| WAV 25.     | Precision alternating current instruments for measuring Volts, Amperes, Watts and Power-factors.   |
| Rapyr 25.   | Radiation pyrometers.  |
| *Bico 25.   | Gerretsen's Bi-colorimeter for Ph. investigations.   |
| *Diaf 25.   | Noyon's Diapherometer.   |
| Mono 25.    | Van Cittert double monochromator.  |
| *Vibra 25.  | Moll vibration galvanometer.   |
| Vacuum 25.  | Moll vacuum thermo-element and vacuum thermo-relais.   |
| *Relais 25. | Moll's thermo-relais.  |


\* Will be published within a few months.

# KIPP-DELFT.

A NEW MOVING COIL GALVANOMETER. Type V.

At present there is a great demand for an inexpensive, simple and robust galvanometer which will entirely answer the expectations of those wanting an instrument which is sufficiently accurate for general laboratory and workshop requirements and suitable for student's use.

The coil, stretched between two strips, carries a plane mirror of 8 mm diameter and moves between the poles of a very powerful steel magnet. Provision is made for clamping the coil and adjusting the zero.

Sensitivity:

1 mm deflection at 1 meter is produced by  $2 \times 10^{-8}$  Amps.  
Resistance about 50 Ohms.  
Full period 3 seconds.

Price only Fl. 90.--.  
( Engl. £ 7.10.0. )  
( U.S.A.\$ 36.--. )

In schools the instrument may conveniently be used for demonstration purposes, e.g. connected to our Moll Small Surface Thermopile, as described in our list "THER 24".

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- |            |  |
|------------|--|
| Galv 27    | Moll Coil Galvanometer.  |
| Ther 24    | Moll Thermopiles.  |
| TS-Z 24    | Torsion String Galvanometer and type Z Galvanometer.                         |
| Ex 24      | Moll Nephelo- and Absorptiometer.  |
| Wav 25     | Precision a.c. instruments for measuring Volts, Amps, Watts and Powerfactor. |
| Reg 25     | Photographically recording drum.   |
| Spectro 25 | Infrared Spectrograph.   |
| Mono 25    | van Cittert Double Monochromator.  |
| Vibra 25   | Moll Vibration Galvanometer.   |
| Relais 25  | Moll Vacuum Thermo-Relay.  |
| Con 27     | Moll Thermo-Converters.  |
| Sol 27     | Solarimeters and Pyrheliometers.   |
| Vacuum 27  | Moll Vacuum Thermoelement with plane window of quartz, etc.                  |
| Tu 29      | Tungsten Arc Lamp Lantern.   |
| Migalv 29  | Moll Micro Galvanometer.   |
| Gezu 29    | Zernike Thermopile.  |
| Vlu 29     | New Instruments of great rapidity for measuring radiation.                   |
| Rooster 28 | Diffraction Gratings.  |
| Mf 29      | Improved Moll Microphotometer.   |

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